ManuZip® volar **Handgelenkbandage**

Wrist Support
Bandage pour poignet
Muñequera
Fasciatura per articolazioni

Art. Nr. 112 220

select

orthopädie
chirurgie
sportmedizin
rehabilitation

TIMBRELLE®
CARE

INHALT/DIRECTORY/SOMMAIRE/ÍNDICE/INDICE

D	Gebrauchsanleitung	4
GB	Directions for use	6
F	Mode d'emploi	8
E	Modo de Empleo	10
I	Istruzione d'uso	12

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Care instructions

Entretien

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Größen

Sizes

Dimensions

Tamaños

Misure 15

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GB Dear Customer,

Thank you for choosing a premium medical product manufactured by BORT GmbH. With the **BORT select® ManuZip® volar Wrist Support**, your specialist medical provider has supplied you with a product that meets the highest medical, functional and quality standards.

All of the products in the BORT **select®** Support series feature Timbrelle antibacterial protection that meets the highest hygienic standards.

The **BORT select® ManuZip® volar Wrist Support** supports the damaged tissue structures. The support's compressive action supports the postural function of the ligaments thus relieving the joint. At the same time, the support helps prevent incorrect movement.

On the move, the compressive knitting and the silicone pads combine to deliver a massage effect, which reduces the tendency for swelling, helps reduce existing swellings and thus accelerates the healing process.

An adjustable volar splint and the **Stabilo®** strap, whose pressure can be adjusted to suit the patient, offer safe stabilization of the wrist.

DIRECTIONS FOR USE

Bend aluminium splint to the desired form before putting it on.

The thumb opening can be opened by unpicking the thread facing the thumb.

The circular **Stabilo®** strap has two differently coloured Velcro surfaces. Always fasten the strap with the transparent area on the fastening point of the brace.

INDICATIONS

For arthrosis of the carpal region, sprains, luxation and functionally conditioned immobilisation. Arthrosis irritation. After removal of plaster cast, arthralgia, drop-hand.

CONTRA INDICATIONS

All medical aids applied externally to the body can, if they are too snugly or tightly fitted, lead to arterial constriction and to oedemas.

Other contraindications include lymph drain malfunctions as well as unexplained soft tissue swelling in areas away from the applied appliance, loss of sensation and circulatory disturbances in the arms and hands as well as skin infections around the treated area of the body.

COMPOSITION

68% Polyamide

25% Elastodies

7% Viscose

IMPORTANT NOTE

No reutilisation – this medical device is meant for the care of one patient only.

The **BORT select® ManuZip® volar Wrist Support** supports the healing process particularly during activities. If instructed to wear the support at night, you must take precautions to avoid impairment of blood circulation. If you experience numbness, please loosen the support immediately and, if necessary, remove it.

If you suffer from persistent complaints or if you experience worsening complaints, you should consult with your doctor.

Please do not apply any creams or ointments at any time while wearing the support as these can destroy the material.

The product contains latex and can trigger allergic reactions.

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